


BOSTON GRAFO ORDER FORM

Date: _____ Due Date: _____ Contact: _____
 Ship To: _____ Account: _____ Phone: _____
 Address: _____ PO#: _____ Fax: _____
 City: _____ State: _____ Zip: _____ Ship Via: _____ E-mail: _____

Patient Name: _____
 Age: _____ Sex: _____ Ht: _____ Wt: _____
 Side: Left Right Bilateral Non-Ambulatory
 Diagnosis: _____

Standard Modifications Are Presumed Unless Otherwise Specified

<p>Ankle Alignment</p> <input type="checkbox"/> Standard- 90° <input type="checkbox"/> As Is <input type="checkbox"/> Correct Ankle to: _____ <input type="checkbox"/> °DF <input type="checkbox"/> °PF <input type="checkbox"/> Feet are asymmetrical	<p>Hindfoot/ Forefoot Alignment</p> <input type="checkbox"/> Standard - Neutral <input type="checkbox"/> Correct FF to: _____ <input type="checkbox"/> °Varus <input type="checkbox"/> °Valgus <input type="checkbox"/> As Is <input type="checkbox"/> Correct HF to: _____ <input type="checkbox"/> °Varus <input type="checkbox"/> °Valgus <input type="checkbox"/> As Is	<p>Add-on Materials</p> <input type="checkbox"/> Molded Inner Boot: 1/8" Opflex (if pad, before pull Standard) <input type="checkbox"/> Heel Post: 1/4" SPP
<p>Solid Ankle GRAFO</p> <input type="checkbox"/> Standard- No Reinforcement <input type="checkbox"/> Foam corregation <input type="checkbox"/> Carbon Reinforcement Specify Other: _____	<p>Padding</p> <input type="checkbox"/> Standard- 1/8" Malleolus Pad + 3/16" Tibial Pad White Colored Foam Before Pull <input type="checkbox"/> Unlined <input type="checkbox"/> After Pull Specify Other: _____	<p>Plastic</p> <input type="checkbox"/> Standard- 3/16" Copoly Specify Other: _____
<p>Trimlines- Encompass Metatarsal STD</p> <input type="checkbox"/> Long Medial <input type="checkbox"/> High Medial <input type="checkbox"/> Long Lateral <input type="checkbox"/> High Lateral Other: _____	<p>Transfer/Color</p> Brace: _____ Straps: _____	<p>Notes</p>



For Avon Use

Mold (s) / Cast (s) Received
 Left Right Bilateral
 Date: _____

Rev 17 4/17