

BOSTON GRAFO ORDER FORM

Date: _____ Due Date: _____ Contact: _____
 Ship To: _____ Account: _____ Phone: _____
 Address: _____ PO#: _____ Fax: _____
 City: _____ State: _____ Zip: _____ Ship Via: _____ E-mail: _____

Patient Name: _____
 Age: _____ Sex: _____ Ht: _____ Wt: _____
 Side: Left Right Bilateral Non-Ambulatory
 Diagnosis: _____

Standard Modifications Are Presumed Unless Otherwise Specified

<p>Ankle Alignment</p> <p><input type="checkbox"/> Standard- 90°</p> <p><input type="checkbox"/> As Is</p> <p><input type="checkbox"/> Correct Ankle to: _____ °DF</p> <p><input type="checkbox"/> _____ °PF</p> <p><input type="checkbox"/> Feet are asymmetrical</p>	<p>Hindfoot/ Forefoot Alignment</p> <p><input type="checkbox"/> Standard - Neutral</p> <p><input type="checkbox"/> Correct FF to: _____ °Varus</p> <p><input type="checkbox"/> _____ °Valgus</p> <p><input type="checkbox"/> As Is</p> <p><input type="checkbox"/> Correct HF to: _____ °Varus</p> <p><input type="checkbox"/> _____ °Valgus</p> <p><input type="checkbox"/> As Is</p>	<p>Add-on Materials</p> <p><input type="checkbox"/> Molded Inner Boot: 1/8" Opflex (if pad, before pull Standard)</p> <p><input type="checkbox"/> Heel Post: 1/4" SPP</p>
<p>Solid Ankle GRAFO</p> <p><input type="checkbox"/> Standard- No Reinforcement</p> <p><input type="checkbox"/> Foam corregation</p> <p><input type="checkbox"/> Carbon Reinforcement</p> <p>Specify Other: _____</p>	<p>Padding</p> <p><input type="checkbox"/> Standard- 1/8" Malleolus Pad + 3/16" Tibial Pad White Colored Foam Before Pull</p> <p><input type="checkbox"/> Unlined <input type="checkbox"/> After Pull</p> <p>Specify Other: _____</p> <p>Trimlines- Encompass Metatarsal STD</p> <p><input type="checkbox"/> Long Medial <input type="checkbox"/> High Medial</p> <p><input type="checkbox"/> Long Lateral <input type="checkbox"/> High Lateral</p> <p>Other: _____</p> <p>Transfer/Color</p> <p>Brace: _____</p> <p>Straps: _____</p>	<p>Plastic</p> <p><input type="checkbox"/> Standard- 3/16" Copoly</p> <p>Specify Other: _____</p> <p>Notes</p>



* Measurements in bold required for finish

For Avon Use

Mold (s) / Cast (s) Received

Left Right Bilateral

Date: _____

Rev 19 8/19