

BOSTON AFO ORDER FORM

Date: _____ Due Date: _____ Contact: _____
 Ship To: _____ Athena # _____ Phone: _____
 Address: _____ PO#: _____ Fax: _____
 City: _____ State: _____ Zip: _____ Ship Via: _____ E-mail: _____

PATIENT INFORMATION

Patient Name: _____ Cast Scan Existing Mold
 Age: _____ Sex: _____ Ht: _____ Wt: _____ Side: Left Right Bilateral Non-Ambulatory

Diagnosis: _____ **Standard Modifications are presumed unless specified**

Ankle Alignment and Correction

Standard 90° L R

Alignment as is L R

Correct Ankle to:

Left °DF _____ °PF _____

Right °DF _____ °PF _____

Heel Post standard dense crepe

Heel Only Extend, taper at toes

Add Heel lift balance to 90 degrees

Match L to R Match R to L

Heel Lift: Left _____ Right _____

Hindfoot/ Forefoot Alignment

**Standard - Neutral Hind foot/
Neutral Forefoot** L R

Hindfoot As Is L R

Left HF Neutral _____ varus _____ valgus

Right HF Neutral _____ varus _____ valgus

Leave Forefoot as is, add post to balance

Left FF Neutral _____ varus _____ valgus

Right FF Neutral _____ varus _____ valgus

Brace Design Additions

Clamshell standard 3/16" aliplast full liner
Provide A/P at calf _____

Tongue: 1/16" LDPE

Wrap Around (3/32" or thinner)

Full Heel Stabilizer 1/4" SPP

Medial heel SPP Lateral Heel SPP

Molded Inner Boot

1/8" Opflex (pad, before pull standard) 3/16" Bock Lite

Padding : Lab Standard- 1/8" Malleolus Pad Matching Colored Foam Before Pull

No Pads in brace Add Pad After Pull

Include Navicular Full Liner 1/8" aliplast

Other _____

Plastic

Lab Standard 1/8" Copoly

Other _____

AFO Design

- Solid
 Semi Solid
 PLS
 Total Contact SAFO

Other _____

Trim Lines: Lab Standard encompass Met Heads

- Proximal Brim Flare
 Full Length Medial High Medial
 Full Length Lateral High Lateral
 Low Profile Trim lines


Transfer/Color/Straps

Brace _____

Straps _____

No color

Other Design Additions



***Measurements in bold required for finish**

AFO Socks:

- Child Small Child Regular
 Adult Small # pairs _____

For Avon Use

Mold (s) / Cast (s) Received

Left Right Bilateral

Date: _____

Rev 2 2019