

# BOSTON AFO ORDER FORM

Date: \_\_\_\_\_ Due Date: \_\_\_\_\_ Contact: \_\_\_\_\_  
 Ship To: \_\_\_\_\_ Athena # \_\_\_\_\_ Phone: \_\_\_\_\_  
 Address: \_\_\_\_\_ PO#: \_\_\_\_\_ Fax: \_\_\_\_\_  
 City: \_\_\_\_\_ State: \_\_\_\_\_ Zip: \_\_\_\_\_ Ship Via: \_\_\_\_\_ E-mail: \_\_\_\_\_

## PATIENT INFORMATION

Patient Name: \_\_\_\_\_  Cast  Scan  Existing Mold  
 Age: \_\_\_\_\_ Sex: \_\_\_\_\_ Ht: \_\_\_\_\_ Wt: \_\_\_\_\_ Side:  Left  Right  Bilateral  Non-Ambulatory

Diagnosis: \_\_\_\_\_ **Standard Modifications are presumed unless specified**

### Ankle Alignment and Correction

**Standard 90°**  L  R

Alignment as is  L  R

Correct Ankle to:

Left  °DF \_\_\_\_\_  °PF \_\_\_\_\_

Right  °DF \_\_\_\_\_  °PF \_\_\_\_\_

### Heel Post standard dense crepe

Heel Only  Extend, taper at toes

Add Heel lift balance to 90 degrees

Match L to R  Match R to L

**Heel Lift:** Left \_\_\_\_\_ Right \_\_\_\_\_

### Hindfoot/ Forefoot Alignment

**Standard - Neutral Hind foot/  
Neutral Forefoot**  L  R

Hindfoot As Is  L  R

Left HF  Neutral \_\_\_\_\_ varus \_\_\_\_\_ valgus

Right HF  Neutral \_\_\_\_\_ varus \_\_\_\_\_ valgus

Leave Forefoot as is, add post to balance

Left FF  Neutral \_\_\_\_\_ varus \_\_\_\_\_ valgus

Right FF  Neutral \_\_\_\_\_ varus \_\_\_\_\_ valgus

### Brace Design Additions

Clamshell standard 3/16" aliplast full liner  
Provide A/P at calf \_\_\_\_\_

Tongue: 1/16" LDPE

Wrap Around (3/32" or thinner)

Full Heel Stabilizer 1/4" SPP

Medial heel SPP  Lateral Heel SPP

### Molded Inner Boot

1/8" Opflex (if pad, before pull standard)  3/16" Bock Lite

### Padding : Lab Standard- 1/8" Malleolus Pad Matching Colored Foam Before Pull

No Pads in brace  Add Pad After Pull

Include Navicular  Full Liner 1/8" aliplast

Other \_\_\_\_\_

### Plastic

Lab Standard 1/8" Copoly

Other \_\_\_\_\_

### AFO Design

- Solid
- Semi Solid
- PLS
- Total Contact SAFO

Other \_\_\_\_\_

### Trim Lines: Lab Standard encompass Met Heads

- Proximal Brim Flare
- Full Length Medial  High Medial
- Full Length Lateral  High Lateral
- Low Profile Trim lines


### Transfer/Color/Straps

Brace \_\_\_\_\_

Straps \_\_\_\_\_

No color

### Other Design Additions



**\*Measurements in bold required for finish**

### AFO Socks:

- Child Small  Child Regular
- Adult Small # pairs \_\_\_\_\_

### For Avon Use

Mold (s) / Cast (s) Received

Left  Right  Bilateral

Date: \_\_\_\_\_

Rev 3 2019