


HOW WILL IT HELP?

STUDIES SHOW BRACING SUCCESS IS RELATED TO WEAR TIME.

The information provided is a report on wear time.

If we are able to reach the prescribed goal, we celebrate success. If we see there are some difficulties, then we can have a conversation and develop strategies to help increase adherence.

This way everyone is working together and no one feels they have to overestimate wear time.



To find a clinic near you, please visit our website.

www.bostonoandp.com

FOLLOW US AND SHARE YOUR STORY!



#BostonO&P



The iButton thermal sensor helps with success



WHAT IS THE IBUTTON?

The iButton is a dime-sized device that periodically takes the inside temperature of the brace to capture a patient's wear patterns.



iButton

As part of our commitment to help patients succeed, we provide the iButton in all Scoliosis braces at no additional charge.




iButton in a Boston Brace

WHY IS WEAR TIME IMPORTANT?

At follow-up visits, the iButton information is downloaded into a simple report.

When patients know they are being monitored, compliance improves.
Extended brace wear improves patient results.

Our goal is to help patients achieve the best possible outcome.



Sample Boston Brace iButton Report

PATIENT/PARENT RESPONSIBILITY

Keeping regularly scheduled follow up appointments helps ensure success. At each follow up, the Boston Brace iButton report will be generated and results reviewed.